VELOPLUS

Wheelchair transport bike



The VeloPlus wheelchair transport bike is designed for transporting people, who can remain seated in their own wheelchair, by cycle. The passenger is not required to be transferred between the wheelchair and the cycle. The co-rider can be wheeled on to the plateau of the cycle easily without having to be lifted. The tilted ramp makes it easy to roll the wheelchair on to the wheel-chair bike and fasten it properly with the wheelchair lock.

The low positioning of the chair on the cycle makes the ride easier and gives a lower centre of gravity. The widely placed and extra steered front-wheels ensure a reduced turning radius and hence, stability. The electric engine may be optionally used for power assist during cycling. The VeloPlus is mostly used in institutions or homes for transporting various persons.

SUITABLE FOR:

The VeloPlus wheelchair bike is suitable for any wheelchair up to 74 cm in width. The wheelchair can be placed quickly and easily by one person on to the cycle and is therefore, suitable for homes or institutions.

TECHNICAL SPECIFICATIONS

Inside leg length Frame height Cycle weight Cycle length Cycle width 80 - 95 cm 52 cm +/- 64 kg 253 cm 110 cm



Hydraulic

disk brakes

X

Wheelchair



/heelchair plateau





